

# forma™

## Allograft Bone



### Allograft Cancellous Crushed

- » 1-4mm and 4-10mm sizes in 5, 10, and 30cc volumes
- » Sterilized with E10<sup>-6</sup> eBeam sterilization technology to preserve native proteins inherent to trabecular tissue
- » Can be used with BMA to create packed bone for osteochondral defects, dental, and lateral gutters in post spinal fusion cases.



### Allograft Demineralized Cube

- » 10x10x10mm, 12x12x12mm, 14x14x14mm sizes
- » Sterilized with E10<sup>-6</sup> eBeam Sterilization Technology
- » Utilization for spinal surgery and dental mandibular canal applications



### Allograft Demineralized Strip

- » Multiple Sizes for use in spinal and dental applications
- » Sterilized with E10<sup>-6</sup> eBeam Sterilization Technology
- » Malleable and flexible to fill lateral gutters and dental and cranio-oral-maxillofacial surgical applications
- » Can be used in conjunction with BMA to make graft osteoconductive



### Allograft Demineralized Paste

- » 1cc, 2.5cc, 5cc and 10cc Syringe
- » Sterilized with R10<sup>-6</sup> Gamma Sterilization Technology
- » 100% Human DBM
- » Easy evacuation and retains a hydrophobic quality during ligation



### Allograft Demineralized Putty

- » 1cc, 2.5cc, 5cc and 10cc Jar
- » Sterilized with R10<sup>-6</sup> Gamma Sterilization Technology
- » 100% DBM with Cortical Cancellous Chips
- » Handling – easy shaped into whatever defect or surgical site repair needed.
- » Hydrophobic qualities during ligation
- » Can be mixed with BMA, PRP to bioactivate proteins to remodel bone

### COMPREHENSIVE SAFETY & QUALITY ASSURANCE PROGRAM

Donor tissue is obtained through a next-of-kin first person consent program with FDA registered recovery agencies contracted with our production partner. These cadaveric donors undergo a series of pre-production serological and microbiological testing, coupled with an extensive eligibility and suitability determination by a medical director prior to release for implant. In addition, the production partner conducts many in-process tests that are congruent with FDA Regulations 21 CFR 1270 and 1271.

**Every donor is screened for potential communicable diseases. These tests include, but are not always limited to the following serological tests:**

- » Hepatitis B core antigen (HBcAg)
- » Hepatitis C antibodies (HCVAb)
- » Human Immunodeficiency Virus 1/O/2 antibodies (HIV-1/O/2 Ab)
- » Human T-lymphotropic virus I/II (HTLVI/II)
- » Anti-Hepatitis B core total antibodies (HBcTotal)
- » RPR – Syphilis Testing

### ADDITIONAL REGULATORY CLASSIFICATIONS

- » All Forma allograft products are processed in ISO Certified Class 5, Class 100 manufacturing environments.
- » All products are tracked in a validated GxP compliant tracking software using serialized inventory controls throughout the entire production through to release for implantation.
- » Forma Allograft is processed by a partner company registered with the United States Food and Drug Administration as a HCT/P establishment and is validated and audited each year.

# forma™

## Allograft Bone

### Crushed Cancellous

PRODUCT NUMBER	DESCRIPTION	SIZE
CS-11015	Cancellous Crush	1-4 mm, 15cc
CS-11030	Cancellous Crush	1-4 mm, 30cc
CS-11060	Cancellous Crush	1-4 mm, 60cc
CS-11115	Cancellous Crush	4-10 mm, 15cc
CS-11130	Cancellous Crush	4-10 mm, 30cc
CS-11160	Cancellous Crush	4-10 mm, 60cc

### Demineralized Cancellous Cube

PRODUCT NUMBER	DESCRIPTION	SIZE
SB-12010	Demineralized Flexible Sponge Cube	10 X 10 X 10 mm
SB-12012	Demineralized Flexible Sponge Cube	12 X 12 X 12 mm
SB-12014	Demineralized Flexible Sponge Cube	14 X 14 X 14 mm

### Demineralized Cancellous Strips

PRODUCT NUMBER	DESCRIPTION	SIZE
SS-12026	Demineralized Flexible Cancellous Strip	26 X 19 X 7 mm
SS-12030	Demineralized Flexible Cancellous Strip	30 X 10 X 7 mm
SS-12031	Demineralized Flexible Cancellous Strip	30 X 20 X 7 mm
SS-12032	Demineralized Flexible Cancellous Strip	30 X 10 X 5 mm
SS-12033	Demineralized Flexible Cancellous Strip	30 X 15 X 5 mm
SS-12050	Demineralized Flexible Cancellous Strip	50 X 20 X 7 mm
SS-12051	Demineralized Flexible Cancellous Strip	50 X 10 X 7 mm
SS-12052	Demineralized Flexible Cancellous Strip	50 X 10 X 5 mm
SS-12053	Demineralized Flexible Cancellous Strip	50 X 15 X 5 mm

### Allograft Bone

PRODUCT NUMBER	DESCRIPTION	SIZE
TB-10001	Patella Bone Block	
TB-10002	Patella Bone Block	Hemi
TB-10003	Iliac Crest Wedge	≤10 mm
TB-10004	Iliac Crest Wedge	11-13 mm
TB-10005	Iliac Crest Wedge	14-17 mm
TB-10006	Iliac Crest Wedge	18-25 mm

### DBM Paste

PRODUCT NUMBER	DESCRIPTION	SIZE
HG-12201	DBM Paste	SYRINGE, 1cc
HG-12202	DBM Paste	SYRINGE, 2.5cc
HG-12205	DBM Paste	SYRINGE, 5cc
HG-12210	DBM Paste	SYRINGE, 10cc

### DBM Putty

PRODUCT NUMBER	DESCRIPTION	SIZE
HG-12301	DBM Putty	VIAL, 1cc
HG-12302	DBM Putty	VIAL, 2.5cc
HG-12305	DBM Putty	VIAL, 5cc
HG-12310	DBM Putty	VIAL, 10cc

### Soft Tissue Allograft

PRODUCT NUMBER	DESCRIPTION	SIZE
ST-14001	Achilles Tendon	W/ Bone Block
ST-14002	Achilles Tendon	W/ Out Bone Block
ST-14003	Achilles Tendon	Pre Sized 10 mm Bone Block
ST-14004	GRACILIS Tendon	30 X 10 X 5 mm
ST-14005	Patella Tendon	Hemi
ST-14006	Patella Tendon	Whole
ST-14007	Patella Tendon	Presized
ST-14008	Peroneus Longus Tendon	
ST-14009	Semitendinosus Tendon	
ST-14010	Tibialis Tendon	Anterior
ST-14011	Tibialis Tendon	Posterior
ST-14013	Fascia Lata	Small
ST-14014	Fascia Lata	Medium
ST-14012	Fascia Lata	Large

MAR-FRM-003

#### ORDERING INFORMATION:

#### Nvision Biomedical Technologies

1350 N Loop 1604 E, Suite 103  
San Antonio, TX 78232

Office: 210-545-3713 ext 107

Email: [Orders@nvisionbiomed.com](mailto:Orders@nvisionbiomed.com)

